REF

OT-1025-8 - Oxy2Mask™ Tyke 7' Tubing
OT-7025-8 - Oxy2Mask™ Tyke Non PVC 7' Tubing

REF

OT-1025-8 - Oxy2Mask™ Tyke 7' Tubing
OT-7025-8 - Oxy2Mask™ Tyke Non PVC 7' Tubing

Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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EC REP logo and contact information: EMERGO EUROPE Westervoortsejdijk 60 6827 AT Arnhem The Netherlands

UK CA 0086 CE 2797 logo and contact information: Emergo Consulting (UK) Limited c/o Cr 360—UL International Compass House, Vision Park Histon Cambridge, CB24 9BZ United Kingdom

Southmedic logo and contact information: 50 Alliance Blvd. Barrie, Ontario Canada L4M 5K3 705-726-9383 705-728-9537 FAX thebetteroxy-

EC REP logo and contact information: EMERGO EUROPE Westervoortsejdijk 60 6827 AT Arnhem The Netherlands

UK CA 0086 CE 2797 logo and contact information: Emergo Consulting (UK) Limited c/o Cr 360—UL International Compass House, Vision Park Histon Cambridge, CB24 9BZ United Kingdom

thebetteroxygenmask.com

* Actual FIO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.
** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

Oxygen Flow Settings	Approximate O ₂ Concentrations
.25 L/min	22%
.50 L/min	25%
.75 L/min	30%
1 L/min	35%
1.5 L/min	43%
2 L/min	47%
3 L/min	50%
≥5 L/min	65%

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