**REF** OKHH-1025-8 - OxyMulti-Mask™ Kid

**OXYGEN MASK WITH HIGH HUMIDITY**

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**REF** OKHH-1025-8 - OxyMulti-Mask™ Kid

**Intended Use:** A minimal contact, open oxygen delivery device for use with USP medical grade oxygen and delivery of aerosolized drug therapy or high humidity therapy.

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**Caution:** Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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**EMERGO EUROPE**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**MedEnvoy Switzerland**  
Gotthardstrasse 28  
6302 Zug  
Switzerland


**UK CA0086 CE 2797**  
**Emergo Consulting (UK) Limited**  
c/o Cr 360—UL International  
Compass House, Vision Park Histon  
Cambridge, CB24 9BZ  
United Kingdom

50 Alliance Blvd.  
Barrie, Ontario  
Canada L4M 5K3  
705-726-9383  
705-728-9537 FAX  
thebetteroxygenmask.com

**EMERGO EUROPE**  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

**MedEnvoy Switzerland**  
Gotthardstrasse 28  
6302 Zug  
Switzerland

**UK CA0086 CE 2797**  
**Emergo Consulting (UK) Limited**  
c/o Cr 360—UL International  
Compass House, Vision Park Histon  
Cambridge, CB24 9BZ  
United Kingdom

Made in Canada

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.


IFU0249J Issue Date: 2024-03-20

Made in Canada

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

IFU0249J Issue Date: 2024-03-20

\*Actual FiO<sub>2</sub> will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.  
\*\* Actual flow rate required should be prescribed in conjunction with pulse oximetry.

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Oxygen Flow Settings	Approximate O <sub>2</sub> Concentrations
1 LPM	22% - 23%
2 LPM	23% - 26%
3 LPM	25% - 28%
4 LPM	28% - 33%
6 LPM	35% - 47%
8 LPM	45% - 58%
10 LPM	53% - 65%
12 LPM	58% - 71%
≥15 LPM	68% - 93%

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