

# OxyMask™ ETCO<sub>2</sub>

Oxygen and End-Tidal CO<sub>2</sub> Therapy



## Patient Safety

- Reduces the risk of CO<sub>2</sub> rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

## Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

## Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

## Operational Efficiencies

One OxyMask™ replaces multiple traditional oxygen therapy interfaces:  
SKU standardization • Inventory reduction • Reduced medical waste  
• Reduces complexity for bedside staff

# OxyMask™ ETCO<sub>2</sub>

- Oxygen flows from 1 to 15+ litres per minute - flush
- Delivers 24% to 65% oxygen with End-Tidal CO<sub>2</sub> monitoring
- Three sizes available:
  - Standard
  - Plus (25% larger mask)
  - Kid (3 - 10 years old, 33 - 70 lbs)
- Does not contain natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OM-2125-8	OxyMask™ End-Tidal CO <sub>2</sub>	Standard	7'	Single
OM-2125-8SLM	OxyMask™ End-Tidal CO <sub>2</sub> , 7' oxygen tubing, 8' sampling line with male connector	Standard	7'	Single
OM-2125-14	OxyMask™ End-Tidal CO <sub>2</sub> with tubing 14'	Standard	14'	Single
OM-2125-14SLM	OxyMask™ End-Tidal CO <sub>2</sub> , 15' tubing, 8' sampling line with male connector	Standard	15'	Single
OC-2025-8	OxyChin™ End-Tidal CO <sub>2</sub> with tubing	Standard	7'	Single
OA-DA-2025	OxyArm™ End-Tidal CO <sub>2</sub> disposable diffuser arm	Standard	None	None
OP-2125-8SLM	OxyPlus™ End-Tidal CO <sub>2</sub> , 7' oxygen tubing, 8' sampling line with male connector	Plus	7'	Single
OP-2325-8	OxyPlus™ End-Tidal CO <sub>2</sub> with tubing	Plus	7'	Single
OP-2125-8	OxyPlus™ End-Tidal CO <sub>2</sub>	Plus	7'	Single
OP-2125-14	OxyPlus™ End-Tidal CO <sub>2</sub> dual elastics & with tubing 14'	Plus	14'	Double
OK-2125-8	OxyKid™ End-Tidal CO <sub>2</sub> with tubing	Kid	7'	Single
OK-2125-8SLM	OxyKid™ End-Tidal CO <sub>2</sub> , 7' oxygen tubing, 8' sampling line with male connector	Kid	7'	Single



Manufactured by:  
 Southmedic Inc.  
 50 Alliance Blvd., Barrie, ON L4M 5K3  
[www.thebetteroxygenmask.com](http://www.thebetteroxygenmask.com)

email: [smedic@southmedic.com](mailto:smedic@southmedic.com)  
 1-800-463-7146 Toll free in North America  
 1-705-726-9383  
 1-705-728-9537 Fax  
 ISO 13485