

OxyMask™

Oxygen Therapy



Patient Safety

- Reduces the risk of CO₂ rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy


Operational Efficiencies

One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
SKU standardization • Inventory reduction • Reduced medical waste

- Reduces complexity for bedside staff

OxyMask™

- Delivers Oxygen flows from 1 to 15+ litres per minute - flush, 24% to 90% FiO₂
(Tyke size 1/4L to 6+, 22% to 65% FiO₂)
- Four sizes available:
 - Standard
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs)
 - Tyke (6 months - 3 years old, 16 - 33 lbs)
- Does not contain natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OM-1125-8	OxyMask™ with 7' tubing	Standard	7'	Single
OM-3025-8	OxyMask™ with Yellow Endotracheal Tube Insert	Standard	7'	Single
OM-1125-14	OxyMask™ with 14' tubing	Standard	14'	Single
OM-1425-8	OxyMask™ with ear elastics and 7' tubing	Standard	7'	Double
OM-7025-8*	OxyMask™ Non PVC with tubing 7'	Standard	7'	Single
OP-1125-8	OxyPlus™ Large Mask w 7' tubing	Plus	7'	Single
OP-3025-8	OxyPlus™ with Yellow Endotracheal Tube Insert	Plus	7'	Single
OP-1125-14	OxyPlus™ Large Mask w 14' tubing	Plus	14'	Single
OP-1425-8	OxyPlus™ Large Mask w ear elastics & 7' tubing	Plus	7'	Double
OP-7025-8*	OxyPlus™ Non PVC with tubing 7'	Plus	7'	Single
OK-1125-8	OxyKid™ Pediatric Mask w 7' tubing	Kid	7'	Single
OK-1125-14	OxyKid™ Pediatric Mask w 14' tubing	Kid	14'	Single
OK-7025-8*	OxyKid™ Non PVC with tubing 7'	Kid	7'	Single
OT-1025-8	OxyTyke™ Pediatric Mask w 7' tubing Non-Painted	Tyke	7'	Single
OT-7025-8*	OxyTyke™ Non PVC with tubing 7'	Tyke	7'	Single

*Non PVC



Manufactured by:
Southmedic Inc.
50 Alliance Blvd., Barrie, ON L4M 5K3
www.southmedic.com
www.thebetteroxygenmask.com

email: smedic@southmedic.com
1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
ISO 13485

01/2016 PO177