

Labour and Delivery

Benefits to the labouring Mom:

Lightweight, open design:

- Reduces feelings of claustrophobia
- Allows unrestricted communication
- Access for fluids and oral medication
- Sits lightly on the face, reducing the possibility of skin breakdown under the nose, on the cheeks and on the ears
- Aids compliance of prescribed oxygen therapy



Benefits for L&D Staff:

Reduces risk of CO₂ rebreathing

- During exhalation, the mask openings allow expired carbon dioxide to escape at flows as low as 1 LPM

Uninterrupted oxygen therapy

- Patient communication, care and oral based therapies are possible without removing the mask and interrupting delivery of oxygen

Reduces risk of contamination and equipment set up error

- By using one OxyMask™ vs. multiple devices, the bedside “clutter” is minimized

Reduces risk of aspiration

- Open mask design minimizes the potential of emesis being aspirated

Time savings

- Time can be better utilized for one-to-one care as oxygen flow can now be quickly titrated to meet the patient’s needs without changing delivery device

Does not contain any natural rubber latex or phthalates such as DEHP

- Minimizes potential for allergic response