

OxyMask™ ETCO₂

Oxygen and End-Tidal CO₂ Therapy



Patient Safety

- Reduces the risk of CO₂ rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

Operational Efficiencies

One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
SKU standardization • Inventory reduction • Reduced medical waste
• Reduces complexity for bedside staff

OxyMask™ ETCO₂

- Oxygen flows from 1 to 12-15 litres per minute
- Delivers 24% to 65% oxygen with End-Tidal CO₂ monitoring
- Three sizes available:
 - Adult
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs)
- Does not contain natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OM-2125-8	OxyMask™ End-Tidal CO ₂	Adult	7'	Single
OM-2325-8	OxyMask™ End-Tidal CO ₂ with tubing with yellow ETTI	Adult	7'	Single
OM-2125-8SLM	OxyMask™ End-Tidal CO ₂ , 7' oxygen tubing, 8' sampling line with male connector	Adult	7'	Single
OM-2125-14	OxyMask™ End-Tidal CO ₂ with tubing 14'	Adult	14'	Single
OM-2125-14SLM	OxyMask™ End-Tidal CO ₂ , 14' tubing, 15' sampling line with male connector	Adult	14'	Single
OC-2025-8	OxyChin™ End-Tidal CO ₂ with 7' oxygen tubing	Adult	7'	Single
OA-DA-2025	OxyArm™ End-Tidal CO ₂ Disposable Diffuser Arm	Adult	None	None
OP-2125-8SLM	OxyMask™ End-Tidal CO ₂ Plus, 7' oxygen tubing, 8' sampling line with male connector	Plus	7'	Single
OP-2325-8	OxyMask™ End-Tidal CO ₂ Plus with tubing with yellow endotracheal tube insert	Plus	7'	Single
OP-2125-8	OxyMask™ End-Tidal CO ₂ Plus with tubing	Plus	7'	Single
OP-2125-14	OxyMask™ End-Tidal CO ₂ Plus dual elastics & with tubing 14'	Plus	14'	Double
OK-2125-8	OxyMask™ End-Tidal CO ₂ Kid with tubing	Kid	7'	Single
OK-2125-8SLM	OxyMask™ End-Tidal CO ₂ Kid, 7' oxygen tubing, 8' sampling line with male connector	Kid	7'	Single



Manufactured by:
 Southmedic Inc.
 50 Alliance Blvd., Barrie, ON L4M 5K3
www.thebetteroxygenmask.com

email: oxymask@southmedic.com
 1-800-463-7146 Toll free in North America
 1-705-726-9383
 1-705-728-9537 Fax
 ISO 13485

02/2017 PO175