

thebetteroxygenmask.com

Adult

Plus
25% larger

Kid
3 - 10 years old,
33 - 70 lbs

Tyke
6 months - 3 years old,
16 - 33 lbs

24% to 90% oxygen concentration - 1 to 15+ lpm to flush

22% to 65% oxygen - ¼ to 5 lpm

OxyMask™
O₂



OxyMask™
OM-1125-8
OM-1125-14
OM-1425-8
OM-7025-8*



OxyMask™ Plus
OP-1125-8
OP-1125-14
OP-1425-8
OP-7025-8*



OxyMask™ Kid
OK-1125-8
OK-1125-14
OK-7025-8*



OxyMask™ Tyke
OT-1025-8
OT-7025-8*

23% to 83% oxygen concentration - 1 to 15+ lpm to flush

22% to 93% oxygen - 1 to 15+ lpm to flush

OxyMulti-Mask™
O₂ / Aerosol



OxyMulti-Mask™
OHH-1025-8
OHH-1425-8
OMMK-4025-8



OxyMulti-Mask™ Plus
OPHH-1025-8
OPHH-1425-8
OPMK-4025-8



OxyUltra™ Kid Plus
OKHH-1025-8
OKHH-1425-8
OKMK-4025-8

23% to 83% oxygen concentration - 1 to 15+ lpm to flush

OxyTrach™
O₂ / Aerosol



OxyTrach™
OAT-1025
OAT-1025-CS

24% to 65% oxygen concentration - 1 to 15 lpm

OxyMask™
ETCO₂
O₂ / ETCO₂



OxyMask™ ETCO₂
OM-2125-8
OM-2325-8
OM-2125-8SLM
OM-2125-14
OM-2125-14SLM



OxyMask™ ETCO₂ Plus
OP-2125-8SLM
OP-2325-8
OP-2125-8
OP-2125-14



OxyMask™ ETCO₂ Kid
OK-2125-8
OK-2125-8SLM

24% to 50% oxygen concentration - 1 to 12-15 lpm

OxyArm™
O₂



OxyArm™
OA-1025-8
OA-DA-1025
OA-DA-3025
OA-DA-1006
OA-DA-3006
OA-DA-2025**
OxyChin™
OC-2025-8**



OxyArm™ Plus
OA-PLUS-1125-8
OA-PLUS-3125-8
OA-PLUS-1125
OA-PLUS-3125
OA-HB-1125



OxyMask™
Aerosol
Aerosol



OxyMask™ Aerosol
OMN-1025



OxyMask™ Aerosol Plus
OPN-1025



OxyMask™ Aerosol Kid
OKN-1025



OxyMask™ Aerosol Tyke
OTN-1025

* Non PVC

** OxyArm / OxyChin with ETCO₂