

OxyMulti-Mask™

High Oxygen • High Humidity • Aerosol Therapy



Patient Safety

- Reduces the risk of CO₂ rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

Operational Efficiencies

- One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
- SKU standardization
 - Inventory reduction
 - Reduced medical waste
 - Reduces complexity for bedside staff

OxyMulti-Mask™

- Oxygen flows from 1 to 15+ litres per minute - flush
- Delivers 23% to 83% FiO₂ (Kid size 22% to 93% FiO₂)
- Three sizes available:
 - Adult
 - Plus (25% larger mask)
 - Kid (3 - 10 years old, 33 - 70 lbs)
- Delivers oxygen through a standard small bore tube, or using a 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O₂ delivery
- Can be used with, or without high humidity or aerosol therapy
- Ideal for acute care/EMS situations when you need high FiO₂, high humidity and/or delivery of medicated aerosol therapy
- Does not contain natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OHH-1025-8	OxyMulti-Mask™ with adapter & 7' tubing	Adult	7'	Single
OHH-1425-8	OxyMulti-Mask™ with ear elastics, adapter & 7' tubing	Adult	7'	Double
OMMK-4025-8	OxyMulti-Mask™ Adult with 22mm adapter & 7' tubing w/nebulizer	Adult	7'	Single
OAT-1025	OxyTrach™ 7' tubing	Adult	7'	Single
OAT-1025-CS	OxyTrach™ 7' tubing with comfort strap	Adult	7'	Single
OPHH-1025-8	OxyMulti-Mask™ Plus with adapter & 7' tubing	Plus	7'	Single
OPHH-1425-8	OxyMulti-Mask™ Plus with ear elastics, adapter & 7' tubing	Plus	7'	Double
OPMK-4025-8	OxyMulti-Mask™ Plus with 22mm adapter & 7' tubing w/nebulizer	Plus	7'	Single
OKHH-1025-8	OxyUltra Kid™ with adapter & 7' tubing	Kid	7'	Single
OKHH-1425-8	OxyUltra Kid™ with ear elastics, adapter & 7' tubing	Kid	7'	Double
OKMK-4025-8	OxyUltra Kid™ with 22mm adapter & 7' tubing w/nebulizer	Kid	7'	Single



Manufactured by:
 Southmedic Inc.
 50 Alliance Blvd., Barrie, ON L4M 5K3
www.thebetteroxygenmask.com

email: oxymask@southmedic.com
 1-800-463-7146 Toll free in North America
 1-705-726-9383
 1-705-728-9537 Fax
 ISO 13485

02/2017 PO176