

OxyTrach™

High Oxygen • High Humidity • Aerosol Therapy



- Oxygen flows from 1 to 15 litres per minute - flush
- Delivers 23% to 83% FiO₂
- Delivers oxygen through a standard small bore tube, or using a 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O₂ delivery
- Can be used with, or without high humidity or aerosol therapy

- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case

Code	Description	Mask Size	Tubing Length	Secure
OAT-1025	OxyTrach with 7' tubing	Adult	7'	Elastic
OAT-1025-CS	OxyTrach with 7' tubing and comfort strap	Adult	7'	Comfort Strap



- Respiratory
- Surgical
- Blades & Scapels
- Dynamic Tissue Systems

Manufactured by:
 Southmedic Inc.
 50 Alliance Blvd.,
 Barrie, ON Canada L4M 5K3
www.thebetteroxygenmask.com

email: oxymask@southmedic.com
 1-800-463-7146 Toll free in North America
 1-705-726-9383
 1-705-728-9537 Fax
 ISO 13485