


# OxyMulti-Mask™

High Oxygen • High Humidity • Aerosol Therapy



## Patient Safety

- Reduces the risk of CO<sub>2</sub> rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

## Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

## Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

## Operational Efficiencies

- One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
- SKU standardization
  - Inventory reduction
  - Reduced medical waste
  - Reduces complexity for bedside staff

# OxyMulti-Mask™

- Oxygen flows from 1 to 15+ litres per minute - flush
- Delivers 23% to 83% FiO<sub>2</sub> (Kid size 22% to 93% FiO<sub>2</sub>)
- Three sizes available:
  - Adult
  - Plus (25% larger mask)
  - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
- Delivers oxygen through a standard small bore tube, or using a 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O<sub>2</sub> delivery
- Can be used with, or without high humidity or aerosol therapy
- Ideal for acute care/EMS situations when you need high FiO<sub>2</sub>, high humidity and/or delivery of medicated aerosol therapy
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OHH-1025-8	OxyMulti-Mask™ with adapter & 7' tubing	Adult	7'	Single
OMMK-4025-8	OxyMulti-Mask™ Adult with 22mm adapter & 7' tubing w/nebulizer	Adult	7'	Single
OAT-1025	OxyTrach™ 7' tubing	Adult	7'	Single
OAT-1025-CS	OxyTrach™ 7' tubing with comfort strap	Adult	7'	Comfort Strap
OPHH-1025-8	OxyMulti-Mask™ Plus with adapter & 7' tubing	Plus	7'	Single
OPMK-4025-8	OxyMulti-Mask™ Plus with 22mm adapter & 7' tubing w/nebulizer	Plus	7'	Single
OKHH-1025-8*	OxyUltra Kid™ with adapter & 7' tubing	Kid	7'	Single

\*OxyUltra Kid™ design not shown



Respiratory  
Surgical  
Blades & Scalpels  
Dynamic Tissue Systems

Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON Canada L4M 5K3  
www.thebetteroxygenmask.com

email: [oxymask@southmedic.com](mailto:oxymask@southmedic.com)  
1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
ISO 13485

P0176 Apr 2018