


# OxyTrach™






High Oxygen  High Humidity  Aerosol Therapy



OxyTrach 

Delivers **23% to 83% FiO<sub>2</sub>** on flows from **1 to 15+** litres per minute



-  Delivers oxygen through a standard small bore tube using the 22mm OD adapter attached to a universal swivel elbow.
-  The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy alongside O<sub>2</sub> delivery
-  Can be used with, or without high humidity or aerosol therapy
-  Elongated design accommodates varying lengths of tracheostomy tubes.
-  OxyTrach™ open design helps minimize the risk of CO<sub>2</sub> rebreathing

## OxyTrach Product Information

Part ID	Description	Packing	Oxygen Tubing Length
OAT - 1025	OxyTrach™ Adult	25 / Box	7'
OAT - 1025 - CS	OxyTrach™ Adult, with Comfort Strap	25 / Box	7'



Explore the full Oxy<sub>2</sub>Mask Family of Products Experience

