
OxyMulti-Mask[™] and OxyTrach[™] Quiz





- 1. What is the FiO2 delivery capability of the OxyMulti-MaskTM? (select one)
 - a) 25% to 50%
 - b) 65% to 90%
 - c) 23% to 83%, (22% to 93% for OxyUltraKid)
- 2. Patient benefits of the OxyMulti-Mask Family include? (select one)
 - a) Reduced feeling of claustrophobia/improved communication
 - b) Ability for oral intake without interrupting oxygen therapy
 - c) Reduction in devices required for oxygen and aerosol therapies
 - d) All of the Above
- 3. There are two OxyMulti-Mask $^{\text{TM}}$ sizes. (select one)
 - a) True
 - b) False
- 4. Which of the following functions can the OxyMulti-MaskTM family of products be used for: (select one)
 - a) The OxyMulti-MaskTM family adapts to a wide variety of aerosol delivery devices
 - b) The OxyMulti-MaskTM can be used to deliver Oxygen and Medicated Aerosol Therapy
 - c) The OxyMulti-MaskTM can be used to deliver High humidity
 - d) All of the above
- 5. Proper placement and fit of the Oxy family of products is not important, as they fit and function like other oxygen masks. (select one)
 - a) True
 - b) False
- 6. What is the proper way to adjust oxygen delivery with the OxyTrachTM? (select one)
 - a) Change the percent setting on the aerosol delivery device
 - b) Adjust the flowmeter setting when using the small bore adapter with oxygen tubing
 - c) All of the above

